

▼ BATTERY CHARGER

BCOB 8A Offboard Battery Charger

▼ BATTERY TRAY AND FRAME SIZE

FS24B Battery Tray for Group 24 Batteries (Total frame width = 25 3/8")

▼ BATTERIES

BAT24PR SigmaTek GEL Batteries – Group 24 (Pair - On chair)

▼ CHOOSE TIRE COLOR (APPLIES TO ALL TIRES OR CASTERS)

GREY

▼ DRIVE WHEEL

- DW14PT 14" Cast Wheel Pneumatic Tire See Dealer
- DW14PI 14" Cast Wheel Pneumatic with Flat-Free Inserts
- DW144PT 4 14" x 4" Extra Wide Cast Wheel Pneumatic Tire See Dealer
- DW144PI 4 14" x 4" Extra Wide Cast Wheel Pneumatic with Flat-Free Inserts See Dealer
- DW12PT 5 12 1/2" Cast Wheel Pneumatic Tire See Dealer
- DW12PI 5 12 1/2" Cast Wheel Pneumatic with Flat-Free Inserts See Dealer

Note 4 The 14" x 4" Extra Wide Wheel increases the overall width of the Base by 2 1/2"
 The option to match the wheels to the shroud colour (FC...M) is not available with the 14" x 4" Extra Wide Wheels.
 The 14" x 4" option is not available with the VR2 non-expandable electronics (ECV902A)

Note 5 12 1/2" Drive wheels reduce overall seat to floor height by 1/4", lowers ground clearance and speed, but increases motor torque

▼ FRONT CASTERS

- FC6S 6" Solid Casters **(Not available with the 14" x 4" Extra Wide Cast Wheel options DW144PT or DW144PI)** See Dealer
- FC8PT 6 8" Pneumatic Casters See Dealer
- FC8PI 6 8" Pneumatic Casters with Flat-Free Inserts

Note 6 The 8" casters are NOT compatible with the 12 1/2" Drive wheel configuration (DW12PT and DW12PI)

▼ SEAT TO FLOOR HEIGHT (Approximate measurement to the top of seat pan ± 1/4")

► Options for a Seat Only

- | | |
|--|---|
| <input type="checkbox"/> STFH155 7 15 1/2" (Possible Seat angle of 0° only) | <input type="checkbox"/> STFH1825 18 1/4" (Possible Seat angle of 0° and 3°) |
| <input type="checkbox"/> STFH1625 16 1/4" (Possible Seat angle of 0° only) | <input type="checkbox"/> STFH1925 19 1/4" (Possible Seat angle of 0°, 3°, 6°) |
| <input type="checkbox"/> STFH1725 17 1/4" (Possible Seat angle of 0° and 3°) | <input type="checkbox"/> STFH2025 20 1/4" (Possible Seat angle of 0°, 3°, 6°, 9°) |

► Options for a Power Tilt or a STD Power Elevating Seat

- | | |
|--|--|
| <input type="checkbox"/> STFH155 7, 8 15 1/2" (Possible Seat angle of 0° only) | <input type="checkbox"/> STFH1775 17 3/4" (Possible Seat angle of 0° only) |
| <input type="checkbox"/> STFH1625 8 16 1/4" (Possible Seat angle of 0° only) | <input type="checkbox"/> STFH1875 18 3/4" (Possible Seat angle of 0° only) |
| <input type="checkbox"/> STFH1725 8 17 1/4" (Possible Seat angle of 0° only) | |

► Options for a Combo - Power Tilt with a STD Power Elevating Seat

- | | |
|---|--|
| <input type="checkbox"/> STFH175 17 1/2" (Possible Seat angle of 0° only) | <input type="checkbox"/> STFH1825 18 1/4" (Possible Seat angle of 0° only) |
|---|--|

Note 7 Requires the 12 1/2" Drive wheels (DW12PT or DW12PI). This lowers ground clearance and speed, but increases motor torque.
Note 8 These seat to floor height measurements are not available with the Heavy-Duty 11" Power Elevating Seat (ADDPESXHD)

▼ SEAT TO FLOOR HEIGHT SPECIFICATION (If not standard)



▼ PRE-SET SEAT FRAME ANGLE

- SFE-0 0° Pre-set seat frame seat angle **(Must be selected when using power options)**
- SFE-3 3° Pre-set seat frame seat angle (Only available with a seat configuration without power options)
- SFE-6 6° Pre-set seat frame seat angle (Only available with a seat configuration without power options)
- SFE-9 9° Pre-set seat frame seat angle (Only available with a seat configuration without power options)

▼ LIGHTING PACKAGE

- EOLP 9 LED Light Package (Includes front beams, side flashers, rear tail light with flashers & brake light) See Dealer

Note 9 Requires the upgraded Joystick with Light Operation or the OMNI (ECVJ6KL, ECELJL, ECRJCL, ECRJCTPL, ECRJC35B or EDIB).

▼ NON-EXPANDABLE ELECTRONICS - PG DRIVES VR2 CONTROLLER AND JOYSTICKS

▼ EXPANDABLE ELECTRONICS - PG DRIVES R-NET CONTROLLERS

- ECEL90 R-net EL 90 Amp Expandable Controller
 - ECR120 R-net PM 120 Amp Expandable Controller See Dealer
 - ECRNH 11 Harness for Expandable Controller
- Note 11** Must be ordered with the expandable controller ECEL90 or ECR120 to interface with Joysticks or OMNI screen

▼ EXPANDABLE ELECTRONICS - PG DRIVES R-NET JOYSTICKS

Selection of one joystick option below is required with the expandable controller ECEL90 or ECR120

- ECO Omit Joystick (For use with an OMNI Display) J
 - ECELJ 12 R-net (LED) Joystick
 - ECRJC 12 R-net Joystick with Color screen J
 - ▶ **JOYSTICKS WITH LIGHTING OPERATION**
 - ECELJL 12 R-net (LED) Joystick with Lighting Operation See Dealer
 - ECRJCL 12 R-net Joystick with Color screen and Lighting Operation See Dealer
 - ▶ **ADVANCED JOYSTICK**
 - ECRJC35B 12 R-net Advanced Bluetooth Joystick (3.5" Screen). Lighting Operation, paddle switches, shortcut keys and stereo jacks See Dealer
- Note 12** This option requires the expandable controller ECEL90 or ECR120

▼ JOYSTICK MOUNT

- JSMF Fixed Joystick Mount
- JSMSA Swing Away Joystick Mount - Height adjustable, Angle adjustable and Multi-Axis See Dealer
- Set-up with standard installation to swing-away sideways Set-up at 35° to swing-away sideways and downwards
- JSMCM Midline Mounted Swing Away Joystick Mount See Dealer

▼ JOYSTICK MOUNT SUPPORT

- JSMFH Joystick mount support Bracket

▼ JOYSTICK HAND MOUNT

- HML Install Joystick on Left hand side HMR Install Joystick on Right hand side

▼ JOYSTICK HANDLE

- JSH Standard Joystick Gimble
- JBAT Stick Handle 3" in Length See Dealer
- JBH Ball Handle See Dealer
- JFBH Foam Ball Handle See Dealer
- JTH T-Handle See Dealer
- JMUSH **New** 2.5" Mushroom Handle See Dealer
- JGP Goalpost Handle See Dealer

▼ ATTENDANT DRIVE CONTROL OPTION

- EADC Attendant Hand Control for R-net EL & PM Controls See Dealer

▼ ATTENDANT JOYSTICK HAND MOUNT

- EADCL Install Attendant Joystick on Rear Left hand side EADCR Install Attendant Joystick on Rear Right hand side See Dealer

▼ ENHANCED DISPLAY & DIGITAL INTERFACE BOX (Included with specialty controls)

- EDIB-F-CA 13 OMNI2 (Enhanced Display & Digital Interface Box with Infrared) J
- Note 13** Allows the switch operation of alternative devices through driver controls. Requires one of the R-net fully Expandable Controllers (ECEL90 or ECR120)

▼ ENHANCED DISPLAY MOUNT

- EDSAM User Swing Away Mount EDML (Install Display on LH side) EDMR (Install Display on RH side) J

▼ ENVIRONMENTAL CONTROL MODULE

<input type="checkbox"/> ECU	14	Environment control module - (Output Module)	See Dealer
Note 14 Requires one of the fully expandable R-net controllers (ECEL90 or ECR120)			

▼ BLUETOOTH MOUSE EMULATOR (Includes USB Dongle)

<input type="checkbox"/> EDBT	15	Bluetooth Mouse Emulator (USB Dongle included)	See Dealer
Note 15 Requires one of the fully expandable R-net controllers (ECEL90 or ECR120) . Not required if the ECRJC35B joystick has been selected Not limited to only one selection. Any two can be installed on the same chair, but only one can be used at a time.			
<input type="checkbox"/> Windows and MAC OS		<input type="checkbox"/> Android Devices	<input type="checkbox"/> iPhone, iPad and iPod

▼ POWER OPTIONS & SEATING

<input type="checkbox"/> SAABA		Amyseat Rehab Seat - Seat Frame only	
<input type="checkbox"/> SAPTXHD	16	Heavy-Duty CG Power Tilt	J
<input type="checkbox"/> SAPRXHD	16, 18	Heavy-Duty Power Recline with Integrated Shear Reduction (ISR)	J
<input type="checkbox"/> SAPTPRXHD	16, 18	Heavy-Duty CG Power Tilt and Power Recline with Integrated Shear Reduction (ISR)	J
<input type="checkbox"/> SAPESXHD	16, 17	Heavy-Duty 11" Power Elevating Seat Module for transfer assist	See Dealer
<input type="checkbox"/> ADDPESXHD	17	Module - Add a Heavy-Duty 11" Power Elevating Seat to either SAPTXHD, SAPRXHD or SAPTPRXHD	See Dealer
Note 16 These options include the fully adjustable Amyseat (SAABA)			
Note 17 The maximum user weight capacity with this option is 400 lbs			
Note 18 It is recommended to use the Contour Back Support option with the power reclines (STDSB or STDSB/CUSTOM)			

▼ ELECTRONICS - SINGLE OR DUAL (1 OR 2 ACTUATORS) POWER SEAT FUNCTION CONTROLLERS

<input type="checkbox"/> EPSEL	19, 20	Tilt or Elevating seat activation through R-NET 90 Amp Expandable Controller (ECEL90)	
<input type="checkbox"/> SAC		Single (1) Power Seat Function Controller to Interface through R-net Electronics	
<input type="checkbox"/> DAC		Dual (2) Power Seat Function Controller to Interface through R-net Electronics	
With the SAC or DAC, please select one activation and position from the list below			
<input type="checkbox"/> SWTJSS	21, 22	Single actuator activation through R-net Joystick or Specialty Controls EPSR2A	
<input type="checkbox"/> SWTJSD	21, 23	Dual actuator activation through R-net Joystick or Specialty Controls EPSRA	
<input type="checkbox"/> SWDUALACT		Dual Activation (Allows to operate all seat functions with external switches and through the joystick)	
<input type="checkbox"/> SWTOGGLE	22 or 23	Single Toggle Switch (Single Mode: DAC - Dual Mode: SAC)	<input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation
<input type="checkbox"/> SWPUSH1	22	Single Push Button Switch (Single Mode: SAC)	<input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation
<input type="checkbox"/> SWPUSH2	22 or 23	Dual Push Button Switch (Single Mode: DAC - Dual Mode: SAC)	<input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation
Note 19 Allows for the activation of a single power option through the VR2 or R-net 90 Amp Controllers. Only the Power Tilt or Elevating seat options are available with this feature.			
Note 20 These options do not require the SAC controller as they connect via the chair's power module directly			
Note 21 Requires the R-net 90 Amp (ECEL90) or the 120 Amp (ECR120) Fully Expandable Controllers			
Note 22 Requires the Single (1 Actuator) Power Seat Function (SAC)			
Note 23 Requires the Dual (2 Actuators) Power Seat Function (DAC)			

▼ ELECTRONICS - MULTIPLE (3 to 6 ACTUATORS) POWER SEAT FUNCTION CONTROLLER

<input type="checkbox"/> MAC		Multiple (3 to 6 Actuators) Power Seat Function Controller to Interface through R-net Electronics	
With the MAC, please select one activation and position from the list below			
<input type="checkbox"/> SWTJSM	24, 25	Multiple Actuator activation through R-net Joystick or Specialty Controls MAC-EPSR	
<input type="checkbox"/> SWDUALACT		Dual Activation (Allows to operate all seat functions with external switches and through the joystick)	
<input type="checkbox"/> SWMULTI	25	4-Way Toggle Switch (Up to 4 Activations)	<input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation
<input type="checkbox"/> SWMULTIT	25	Triple Toggle Switch (Up to 5 Activations)	<input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation
<input type="checkbox"/> SWPUSH4	25	4-Way Push Button Switch	<input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation
<input type="checkbox"/> SWPUSH5	25	5-Way Push Button Switch	<input type="checkbox"/> User Activation <input type="checkbox"/> Attendant Activation
Note 24 Requires the Multiple (3 to 6 Actuators) Power Seat Function Controller (MAC)			
Note 25 Requires the R-net 90 Amp (ECEL90) or 120 Amp (ECR120) Fully Expandable Controller, and the multiple Actuator Controller (MAC)			

▼ SWITCH ACTIVATION HAND MOUNT

Switch Installation: TSML (Install Switch on LH side) TSMR (Install Switch on RH side)

▼ LARGE ADJUSTABLE SEAT WIDTH RANGE BETWEEN 19" AND 26"

- SOWN19-24-22 Width adjustment from 19" to 24" with a maximum Depth 22"
- SOWN24-26-22 Width adjustment from 24" to 26" with a maximum Depth 22"

- | | | |
|--|---------------------------------------|--|
| <input type="checkbox"/> 19" (SOWN19L) | <input type="checkbox"/> 22" (SOWN22) | <input type="checkbox"/> 24" (SOWN24HDL) |
| <input type="checkbox"/> 20" (SOWN20L) | <input type="checkbox"/> 23" (SOWN23) | <input type="checkbox"/> 25" (SOWN25HDL) |
| <input type="checkbox"/> 21" (SOWN21) | <input type="checkbox"/> 24" (SOWN24) | <input type="checkbox"/> 26" (SOWN26HDL) |

▼ ADJUSTABLE SEAT DEPTH RANGE BETWEEN 17" AND 22"

- | | |
|---------------------------------------|---------------------------------------|
| <input type="checkbox"/> 17" (SODN17) | <input type="checkbox"/> 18" (SODN18) |
| <input type="checkbox"/> 19" (SODN19) | <input type="checkbox"/> 20" (SODN20) |
| <input type="checkbox"/> 21" (SODN21) | <input type="checkbox"/> 22" (SODN22) |

▼ BACK CANE HEIGHT (The back cane height measurement is taken from the seat pan to the top of the push handles)

- | | |
|--|---|
| <input type="checkbox"/> 16" (SOBCH16) | <input type="checkbox"/> 22" (SOBCH22) 26 |
| <input type="checkbox"/> 18" (SOBCH18) | <input type="checkbox"/> 24" (SOBCH24) |
| <input type="checkbox"/> 20" (SOBCH20) | <input type="checkbox"/> 26" (SOBCH26) |

Note 26 22" is the minimum back cane height required with a power or manual recline.

▼ PRE-SET SEAT ANGLE (For Non-Reclining Backs)

- SOBAA Please Select one: 86° 90° 94° 98° 102° 106° 110° 114° 118° 122°

▼ BACK SUPPORT

- | | | |
|---------------------------------------|---|------------|
| <input type="checkbox"/> STDSB | Standard contour back support 19" to 20" Wide with a maximum shell height of 22" | See Dealer |
| <input type="checkbox"/> STDSB/CUSTOM | Standard contour back support Non Standard Sizes up to 26" wide Width: _____ Depth: _____ | See Dealer |
| <input type="checkbox"/> | Back Upholstery | |
| <input type="checkbox"/> | Omit Back Upholstery | |

Note 27 It is recommended to use the Back Support option with the power recline.

▼ BACK SUPPORT TOTAL HEIGHT (MEASURED FROM THE SEAT PAN TO THE TOP OF THE BACK)

- | | | |
|---------------------------------------|--|---------------------------------------|
| <input type="checkbox"/> 16" (SOBH16) | <input type="checkbox"/> 20" (SOBH20) | <input type="checkbox"/> 24" (SOBH24) |
| <input type="checkbox"/> 18" (SOBH18) | <input type="checkbox"/> 22" (SOBH22) 28 | <input type="checkbox"/> 26" (SOBH26) |
| | | <input type="checkbox"/> 28" (SOBH28) |

Note 28 22" is the minimum back height required with a power or manual recline.

Note 29 The space between the bottom of the back and the seat pan is 4" with a non-reclining back and 6" with a reclining back

▼ SEAT CUSHIONS

▼ CANTILEVER - HEIGHT ADJUSTABLE ARMRESTS - OPENS WITH A 5° ANGLE (For Non-Reclining Backs - Choose 2)

- | | | |
|-----------------------------------|---|---|
| <input type="checkbox"/> SOAFBFLH | Flip-Back Cantilever - Full Length Left Side (Height Adjustable from 7" to 15" with 1" increments) | J |
| <input type="checkbox"/> SOAFBFRH | Flip-Back Cantilever - Full Length Right Side (Height Adjustable from 7" to 15" with 1" increments) | J |
| <input type="checkbox"/> SOAFBDLH | Flip-Back Cantilever - Desk Length Left Side (Height Adjustable from 7" to 15" with 1" increments) | J |
| <input type="checkbox"/> SOAFBDRH | Flip-Back Cantilever - Desk Length Right Side (Height Adjustable from 7" to 15" with 1" increments) | J |

▼ DUAL POST - HEIGHT ADJUSTABLE ARMRESTS - INCLUDES TRANSFER GRIPS (For Non-Reclining Backs - Choose 2)

- | | |
|-----------------------------------|---|
| <input type="checkbox"/> SOADPFLH | Dual Post Height Adjustable - Full Length Left Side (Height Adjustable from 10 ¼" - 15 ¼") |
| <input type="checkbox"/> SOADPFRH | Dual Post Height Adjustable - Full Length Right Side (Height Adjustable from 10 ¼" - 15 ¼") |
| <input type="checkbox"/> SOADPDLH | Dual Post Height Adjustable - Desk Length Left Side (Height Adjustable from 10 ¼" - 15 ¼") |
| <input type="checkbox"/> SOADPDRH | Dual Post Height Adjustable - Desk Length Right Side (Height Adjustable from 10 ¼" - 15 ¼") |

▼ RECLINING FLIP-UP, HEIGHT ADJUSTABLE ARMRESTS (For Reclining Backs - Choose 2)

- | | | |
|-----------------------------------|---|----|
| <input type="checkbox"/> SOANIFLH | Reclining armrest - Full Length Left Side (Height Adjustable from 10 ½" to 14" with ½" increments) | J* |
| <input type="checkbox"/> SOANIFRH | Reclining armrest - Full Length Right Side (Height Adjustable from 10 ½" to 14" with ½" increments) | J* |
| <input type="checkbox"/> SOANIDLH | Reclining armrest - Desk Length Left Side (Height Adjustable from 10 ½" to 14" with ½" increments) | J* |
| <input type="checkbox"/> SOANIDRH | Reclining armrest - Desk Length Right Side (Height Adjustable from 10 ½" to 14" with ½" increments) | J* |
- *Only available with approval for power recline (MWP)

▼ ARMREST HEIGHT ADJUSTMENT AT FACTORY

Please set the Armrest Height at: _____ " (The Armrest Height is taken from the seat pan to the top of the armpad. Not required for Dual Post)

▼ ARMPADS

► Foam Type Selection

<input type="checkbox"/>	SOFAPD	Desk Length Standard (10")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	
<input type="checkbox"/>	SOFAPF	Full Length Standard (14")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	
<input type="checkbox"/>	SOWAPD	30 Desk Length Waterfall (10")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	
<input type="checkbox"/>	SOWAPF	30 Full Length Waterfall (14")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	

► Hybrid Gel and Foam Type Selection (The Gel always stays soft)

<input type="checkbox"/>	SOGELD	Desk Length Hybrid Gel and Foam Flat armpad (10" x 2.5")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	See Dealer
<input type="checkbox"/>	SOGELF	Full Length Hybrid Gel and Foam Flat armpad (14" x 2.5")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	See Dealer
<input type="checkbox"/>	SOWGELD	30 Desk Length Hybrid Gel and Foam Waterfall armpad (10")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	See Dealer
<input type="checkbox"/>	SOWGELF	30 Full Length Hybrid Gel and Foam Waterfall armpad (14")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	See Dealer

► Specialty Type Selection

<input type="checkbox"/>	SOCAR	31 Arm Trough with Hand Support (18" x 4 1/2")	Select Side: <input type="checkbox"/> Left Hand Side	<input type="checkbox"/> Right Hand Side	See Dealer
--------------------------	-------	--	--	--	------------

Note 30 The Waterfall armpads may interfere with flipping back the armrests when combined with Lateral Trunk supports
 Note 31 Arm Trough Armpads may interfere with Armrests Flipping back.

▼ FOOTRESTS

<input type="checkbox"/>	SOAAF	35 Heavy Duty Swing-Away footrests	Set at: <input type="checkbox"/> 60°	<input type="checkbox"/> 70°	
<input type="checkbox"/>	ADDPLHD	32, 35 Heavy Duty Power Elevating & Articulating Legrests	Select Footrest Range: <input type="checkbox"/> SOFLHD12 (12 1/2" - 15 1/2")	<input type="checkbox"/> SOFLHD14 (14 1/2" - 19 1/2")	J*
		35 Please set the Footrest Length at: _____ "	Please Select Function: <input type="checkbox"/> Individual ³³	<input type="checkbox"/> Combined ³⁴	
			Select Footrest Range: <input type="checkbox"/> SOFLPW12 (12 1/2" - 15 1/2")	<input type="checkbox"/> SOFLPW14 (14 1/2" - 19 1/2")	

Note 32 This option includes calf pads
 Note 33 Individual legs count as two power functions.
 Note 34 Combined legs count as one power function.
 Note 35 The Footrest Length is taken from the top of the seat pan to the top of the footplate.

▼ FOOTPLATES & ACCESSORIES

<input type="checkbox"/>	SOSFAF	Flip-up Footplates			
<input type="checkbox"/>	SOAFP	Angle Adjustable & Flip-up Footplates - Each			J
<input type="checkbox"/>	SOAHL	Heel Loops - Each			
<input type="checkbox"/>	SOACPFRR	Calf Pads for Adult Heavy-Duty Swing-Away footrests (SOAAF) with adjustable hardware - Each			See Dealer

▼ CENTER MOUNTS

<input type="checkbox"/>	SOAAF90	38, 39, 40 90° Fixed Center Mount Interface Bracket (Offers Manual adjustability from 95° to 70°)	Select Foot Platform Range: <input type="checkbox"/> SOFBA12 (12 1/2" - 16 1/2")	<input type="checkbox"/> SOFBA16 (16 1/2" - 19 1/2")	
<input type="checkbox"/>	ADDPFHD	36, 37, 38, 39, 40 90° Power Center Mount Elevating & Articulating Foot Platform	Select Foot Platform Range: <input type="checkbox"/> SOFBPW12 (12 1/2" - 16 1/2")	<input type="checkbox"/> SOFBPW16 (16 1/2" - 20 1/2")	J*

Please set the Foot Platform Length at: _____ "
 *Only available with approval for power ELR's (MWP)
 Note 36 The power center mount counts as one power function
 Note 37 This option includes calf pads
 Note 38 Links and attaches the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB) to the seat
 Note 39 Please select the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB)
 Note 40 The Foot Platform length is taken from the top of the seat pan to the top of the foot platform or split footplates.

▼ FOOT PLATFORMS & ACCESSORIES

<input type="checkbox"/>	SOFUFP	Flip-up one piece Foot Platform			J
<input type="checkbox"/>	SOAFPFB	Angle Adjustable & Flip-up Split Footplates for center mount - Each			J
<input type="checkbox"/>	SOACPCM	Calf Pads for 90° Fixed Center Mount (SOAAF90) with adjustable hardware - Each			See Dealer

▼ SEAT POSITIONING (NOT COMPATIBLE WITH THE WC19 ----- ISO OPTION)

- SOAUTO80 80" Seat Positioning Strap

▼ HEADREST

- Whitmyer P1L Linx Headrest
 M2300 Bracket
 SOHRSB Bracket to Mount on the Seat Frame crossbar (**Without Recline**)
 SOHRPR Bracket to Mount on the Seat Frame (**For Reclining Backs**)

Please select one of the headrest options with the mounting hardware and specify where the headrest will be installed (The frame or back cushion)

▼ LATERAL TRUNK SUPPORTS

- SOSLS3-C Adjustable Laterals with Small Curved Pads (Each) - 3" x 5" See Dealer
 SOSLM4-C Adjustable Laterals with Medium Curved Pads (Each) - 4" x 6" See Dealer
 SOSLL5-C Adjustable Laterals with Large Curved Pads (Each) - 5" x 7" See Dealer

▼ HIP & THIGH SUPPORTS

- SOHTSS10 Adjustable Hip & Thigh Support Pads (Each) - Small 10" x 4" See Dealer
 SOHTSM12 Adjustable Hip & Thigh Support Pads (Each) - Medium 12" x 5" See Dealer
 SOHTSL14 Adjustable Hip & Thigh Support Pads (Each) - Large 14" x 6" See Dealer
 SOHTSARH Adjustable and Removable Hip & Thigh Support Hardware - Each See Dealer

Please select one of the hip & thigh support pads and the adjustable & removable hardware.

▼ ACCESSORIES

- SOCH Cup Holder (Self-Leveling) See Dealer
 SOCPH Ram Mounts X-Grip Cell Phone Holder See Dealer
 SOTHS Ram Mounts X-Grip Small Tablet Holder See Dealer
 SOTHL Ram Mounts X-Grip Large Tablet Holder with Swivel See Dealer

▼ REFLECTOR KIT

- REFLECT Red & White Reflector Kit Stickers (They will not be installed - Included in a bag with the owner's manual) See Dealer

▼ SWITCHES AND OPTIONS

▼ SPECIALTY DRIVER CONTROL OPTIONS

▼ Therafin Sip n' Puff Kit

<input type="checkbox"/>	SCSNPK	41,42	Therafin Sip n' Puff Kit			J
<input type="checkbox"/>	SCSNPBTk	42	Therafin Sip n' Puff Breath Tube Kit			J
<input type="checkbox"/>	SCSRMSNP	42	Stealth Rock-n-Lock Swingaway Arm for Sip-N-Puff	<input type="checkbox"/> SWL (Swings on Left)	<input type="checkbox"/> SWR (Swings on Right)	J
<input type="checkbox"/>	SCSRM	42	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings on Left)	<input type="checkbox"/> SWR (Swings on Right)	J
<input type="checkbox"/>	SSES-STD	42, 44	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	J
<input type="checkbox"/>	SSBS-STD	42	Buddy Button Switch (Stop Switch) - Red			J
<input type="checkbox"/>	SOMAH	42	Multi-axis Headrest Removable Mounting Hardware			J
<input type="checkbox"/>	SOFHRS	43	Small Curved Headrest Pad 4" x 9"			J
<input type="checkbox"/>	SOFHRM	43	Medium Curved Headrest Pad 5" x 10"			J
<input type="checkbox"/>	SOFHRL	43	Large Curved Headrest Pad 6" x 11"			J

▼ Switch-It Proximity Head Array

<input type="checkbox"/>	SCSIHAA	41, 42	Adult - 3 Switch Proximity Head Array with 1/8" Mono Jack for mode/reverse			See Dealer
<input type="checkbox"/>	SCSRM	42	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings on Left)	<input type="checkbox"/> SWR (Swings on Right)	See Dealer
<input type="checkbox"/>	SSES-STD	42, 44	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	See Dealer
<input type="checkbox"/>	SSBS-STD	42	Buddy Button Switch (Stop Switch) - Red			See Dealer
<input type="checkbox"/>	SOMAH	42	Multi-axis Headrest Removable Mounting Hardware			See Dealer

▼ Switch-It Proximity Head Array & Sip-N-Puff Combo

<input type="checkbox"/>	SCSICMBOA	41,42	Adult-4 Switch Combo Head Array with two proximity sensors for Left & Right and Low pressure Sip-N-Puff for forward and reverse and 1/8" jack for			See Dealer
<input type="checkbox"/>	SCSNPBTk	42	Therafin Sip-N-Puff Breath Tube Kit			See Dealer
<input type="checkbox"/>	SCSRMSNP	42	Stealth Rock-n-Lock Swingaway Arm for Sip-N-Puff	<input type="checkbox"/> SWL (Swings on Left)	<input type="checkbox"/> SWR (Swings on Right)	See Dealer
<input type="checkbox"/>	SCSRM	42	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings on Left)	<input type="checkbox"/> SWR (Swings on Right)	See Dealer
<input type="checkbox"/>	SSES-STD	42, 44	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	See Dealer
<input type="checkbox"/>	SSBS-STD	42	Buddy Button Switch (Stop Switch) - Red			See Dealer
<input type="checkbox"/>	SOMAH	42	Multi-axis Headrest Removable Mounting Hardware			See Dealer

▼ Proportional RIM Head Control Package

<input type="checkbox"/>	SCPHC	41, 42	Proportional RIM Head Control Package			J
<input type="checkbox"/>	SCSRM	42	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings on Left)	<input type="checkbox"/> SWR (Swings on Right)	J
<input type="checkbox"/>	SSES-STD	42,44	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	J
<input type="checkbox"/>	SOMAH	42	Multi-axis Headrest Removable Mounting Hardware			J

▼ Switch-It MicroPilot Light Touch Proportional Mini joystick for chin / lip control

<input type="checkbox"/>	SCSIMP	41, 42	MicroPilot Light touch mini Proportional Joystick			See Dealer
<input type="checkbox"/>	SCSRMMP	42	Stealth Rock-n-Lock Swingaway Arm for Chin MP	<input type="checkbox"/> SWL (Swings on Left)	<input type="checkbox"/> SWR (Swings on Right)	See Dealer
<input type="checkbox"/>	SCSRM	42	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	<input type="checkbox"/> SWL (Swings on Left)	<input type="checkbox"/> SWR (Swings on Right)	See Dealer
<input type="checkbox"/>	SSES-STD	42, 44	Egg Switch for Changing Modes - Black	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	See Dealer
<input type="checkbox"/>	JCC		Chin Cup			See Dealer
<input type="checkbox"/>	SOMAH	42	Multi-axis Headrest Removable Mounting Hardware			See Dealer
<input type="checkbox"/>	SOFHRS	43	Small Curved Headrest Pad 4" x 9"			See Dealer
<input type="checkbox"/>	SOFHRM	43	Medium Curved Headrest Pad 5" x 10"			See Dealer
<input type="checkbox"/>	SOFHRL	43	Large Curved Headrest Pad 6" x 11"			See Dealer

▼ Proportional Compact or Mini joystick for Midline Mount

<input type="checkbox"/>	SCPGMJS	41, 42	PG Drives JC200 Proportional Compact Remote Joystick for Midline Mount			See Dealer
<input type="checkbox"/>	SCGMJSM	45	Gatlin Midline Swing-away Mount for JC200	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	See Dealer
<input type="checkbox"/>	SCSIMP	41, 42	Switch-It MicroPilot Light touch mini Proportional Joystick			See Dealer
<input type="checkbox"/>	SCGMPM	46	Gatlin Midline Swing-away Mount for MicroPilot	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	See Dealer
<input type="checkbox"/>	SCBTM	47	Bullet Hand Tray for MicroPilot	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	See Dealer

▼ Proportional Compact or Mini joystick for Armrest Mount

<input type="checkbox"/>	SCPGMJS	41, 42	PG Drives JC200 Proportional Compact Remote Joystick for Armrest Mount			See Dealer
<input type="checkbox"/>	SCARMJSM	45	Swing-away and Multi-Axis armrest mount	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	See Dealer
<input type="checkbox"/>	SCSIMP	41, 42	Switch-It MicroPilot Light touch mini Proportional Joystick			See Dealer
<input type="checkbox"/>	SCARMPM	46	Swing-away and Multi-Axis armrest mount	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	See Dealer
<input type="checkbox"/>	SCBTM	47	Bullet Hand Tray for MicroPilot	<input type="checkbox"/> MLH (Mount on Left)	<input type="checkbox"/> MRH (Mount on Right)	See Dealer

Note 41 Requires the OMNI Enhanced Display and Digital Interface Box (EDIB)

Note 42 This option must be selected for the kit

Note 43 Choose one (1) Headrest Pad

Note 44 The mode switch type can be changed from this selection on Page #6

Note 45 Must be selected for use with the PG Drives JC200 Mini Joystick

Note 46 Must be selected for use with the Switch-It MicroPilot

Note 47 This option is only compatible with the Switch-It MicroPilot

